



Getein
Biotech, Inc.

Stock Code: 603387

OPTIMIZED POINT-OF-CARE SOLUTION

MAKE TEST EASY

Getein 1100

Immunofluorescence Quantitative Analyzer



Getein 1100

Immunofluorescence Quantitative Analyzer



HIGHLY EFFICIENT & ACCURATE

Advanced fluorescence immunoassay

Multiple quality control



REAL-TIME AND RAPID TEST


One-step test

3-15 min/test


5 sec/test for multiple tests

OPERATION MODES

Inside Mode (single sample rapid test mode)



Sample Dispense



Test Card Insert



Click "Start" Icon




Result Show and Print

Quick mode (mass samples rapid test mode)



Sample Dispense



Timing the Reaction Manually



Click "Start" Icon



Result Show and Print



CONVENIENT OPERATION

RFID card calibration

Keyboard and mouse connectivity through USB port

Handwriting input available

Continuous test for 3 hours with lithium battery



USER-FRIENDLY INTERFACE

Android system

7 inch touch screen



TECHNICAL PARAMETERS

Methodology

Lateral flow chromatography
(Immunofluorescence)

Result

Quantitative

Sample Type

WB, plasma, serum, urine

Storage Capacity

10000 data

Language

English/Chinese/Spanish/Portuguese

Screen

7 inch touch screen

Power Supply

100-240V~50Hz/60Hz, 60VA

Working Environment

Temperature: 10-35 °C
Relative humidity ≤70%
Air pressure 70.0~106.0kpa

Dimension

261mm×241mm×115mm (D×W×H)

Weight

2.0 kg

TEST ITEMS

Cat. #	TEST ITEMS	DISEASES	CUT-OFF VALUE	SAMPLE TYPES	MEASURING RANGE	SAMPLE VOLUME	REACTION TIME	QUALIFICATION
Cardiac Markers								
IF1001	cTnI	Myocardial infarction	0.1ng/ml	S/P/WB	0.1-50.0ng/ml	100µL	10min	CFDA CE
IF1002	NT-proBNP	Heart failure	300pg/ml	S/P/WB	100-35000pg/ml	100µL	10min	CFDA CE
IF1004	NT-proBNP/cTnI	Heart failure/ acute coronary syndrome	NT-proBNP: 300pg/ml cTnI: 0.1ng/ml	S/P/WB	100-15000pg/ml 0.1-50.0ng/ml	100µL	10min	CFDA CE
IF1005	CK-MB/cTnI/Myo	Myocardial damage /infarction	CK-MB: 5.0ng/ml cTnI: 0.1ng/ml Myo: 70ng/ml	S/P/WB	2.5-80.0ng/ml 0.1-50.0ng/ml 30.0-600.0ng/ml	100µL	10min	CFDA CE
IF1006	D-Dimer	Venous thromboembolism	0.5mg/L	P/WB	0.1-10.0mg/L	100µL	10min	CFDA CE
IF1012	CK-MB/cTnI	Myocardial damage /infarction	CK-MB: 5.0ng/ml cTnI: 0.1ng/ml	S/P/WB	2.5-80.0ng/ml 0.1-50.0ng/ml	100µL	10min	CE
IF1014	H-FABP	Myocardial damage	6.36ng/ml	S/P/WB	1.0-120.0ng/ml	100µL	3min	CFDA CE
IF1016	CK-MB/cTnI/H-FABP	Myocardial damage /infarction	CK-MB: 5.0ng/ml cTnI: 0.1ng/ml H-FABP: 6.36ng/ml	S/P/WB	2.5-80.0ng/ml 0.1-50.0ng/ml 2.0-100.0ng/ml	100µL	10min	CFDA CE
IF1018	CK-MB	Myocardial injury	5.0ng/ml	S/P/WB	2.5-80.0ng/ml	100µL	10min	CE
Diabetes Mellitus								
IF1017	HbA1c	Diabetes mellitus	3.8%-5.8%	WB	2%-14%	10µL	5min	NGSP CFDA IFCC CE
Inflammation								
IF1003	hs-CRP+CRP	Cardiovascular inflammation /normal inflammation	3mg/L 10mg/L	S/P/WB/ Fingertip Blood	0.5-200mg/L	10µL	3min	CFDA CE
IF1007	PCT	Sepsis, bacterial infection	0.1ng/ml	S/P/WB	0.1-50.0ng/ml	100µL	15min	CFDA CE
IF1015	PCT/CRP	Sepsis, bacterial infection	PCT: 0.1ng/ml CRP: 3.0mg/L	S/P/WB	0.1-50.0ng/ml 0.5-200.0mg/L	100µL	15min	CFDA CE
IF1044	SAA	Bacterial/Virus infection	10.0mg/L	S/P/WB/ Fingertip Blood	5.0-200.0mg/L	10µL	5min	CE
Renal Function								
IF1008	CysC	Acute and chronic renal diseases	0.51-1.09mg/L	S/P/WB	0.5-10.0mg/L	10µL	3min	CFDA CE
IF1009	mAlb	Diabetic nephropathy, hypertensive nephropathy	20.0mg/L	Urine	10.0-200.0mg/L	100µL	3min	CFDA CE
IF1010	NGAL	Acute kidney injury	Serum: 200ng/ml Urine: 100ng/ml	S/Urine	50-5000ng/ml	10µL	10min	CFDA CE
IF1011	β ₂ -MG	Acute and chronic kidney diseases/tumours	0.8-3.0mg/L	S/P/WB	0.5-20.0mg/L	10µL	3min	CFDA CE
Hormone								
IF1024	TSH	Thyroid malfunction	0.27-4.20µIU/ml	S/P	0.10-50.00µIU/ml	100µL	15min	CE
IF1022	T3	Hyperthyroidism, hypothyroidism	1.30-3.10nmol/L	S/P	0.30-10.00nmol/L	40µL	15min	CE
IF1023	T4	Hyperthyroidism, hypothyroidism	59.0-154.0nmol/L	S/P	5.4-320.0nmol/L	40µL	15min	CE
Tumor								
IF1053	tPSA	Prostate cancer	4.0ng/mL	S/P	0.50ng/mL-100.00ng/mL	100µL	15min	
Bone Metabolism								
IF1031	25-OH-VD	Osteomalacia, osteoporosis	30.0-50.0ng/mL	S/P	8.0-70.0ng/mL	40µL	15min	CE
Fertility								
IF1013	HCG+β	Fertility	5.1mIU/ml	S/P	5-100000mIU/ml	100µL	10min	CFDA CE

Coming Soon: FOB, LH, FSH, PRL, AMH...

GP Getein Biotech, Inc.

Add: No.9 Bofu Road, Luhe District, Nanjing, 211505, China
Tel: +86-25-68568508/68568594
Fax: +86-25-68568500
E-mail: sales@getein.com.cn; overseas@getein.com.cn
Web: www.bio-GP.com.cn

FSC CE CFDA NGSP IFCC

